

KERN 357-07

KERN

Individual weight, OIML Class M2, 100 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-637

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	357

Measuring System

Test weight - nominal value	100 g
Test weight - volume (cm ³)	11,9048
Test weight - volume uncertainty (cm ³)	0,2409
Permissible tolerance according to OIML (+/-)	16 mg
Test weight - Density [Min] (OIML)	2300
Weight OIML class	M2

Construction

Shape weight	Knob
Dimensions (Ø×H)	22×37 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	24×24×51 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg