

KERN 357-08

KERN

Individual weight, OIML Class M2, 200 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-638

Pictograms

STANDARD



Category

| | |
|------------------|-------------------|
| Brand | KERN |
| Product category | Individual weight |
| Product group | OIML M2 |
| Product family | 357 |

Measuring System

| | |
|---|---------|
| Test weight - nominal value | 200 g |
| Test weight - volume (cm ³) | 23,8095 |
| Test weight - volume uncertainty (cm ³) | 0,4819 |
| Permissible tolerance according to OIML (+/-) | 30 mg |
| Test weight - Density [Min] (OIML) | 2300 |
| Weight OIML class | M2 |

Construction

| | |
|------------------|-----------------------------------|
| Shape weight | Knob |
| Dimensions (Ø×H) | 28×47 mm |
| Material | stainless steelfine turned (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 31×31×56 mm |
| Net weight | 0,2 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,20 kg |
| Gross weight approx. | 0,20 kg |
| Shipping weight | 0,2 kg |