

# KERN 357-09

**KERN**

Individual weight, OIML Class M2, 500 g, Knob, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration 962-639

## Pictograms

## STANDARD



## Category

|                  |                   |
|------------------|-------------------|
| Brand            | KERN              |
| Product category | Individual weight |
| Product group    | OIML M2           |
| Product family   | 357               |

## Measuring System

|   |         |
|---|---------|
| Test weight - nominal value                         | 500 g   |
| Test weight - volume (cm <sup>3</sup> )             | 59,5238 |
| Test weight - volume uncertainty (cm <sup>3</sup> ) | 1,2046  |
| Permissible tolerance according to OIML (+/-)       | 80 mg   |
| Test weight - Density [Min] (OIML)                  | 2300    |
| Weight OIML class                                   | M2      |

## Construction

|                  |                                   |
|------------------|-----------------------------------|
| Shape weight     | Knob                              |
| Dimensions (Ø×H) | 38×63 mm                          |
| Material         | stainless steelfine turned (OIML) |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 57×57×74 mm    |
| Net weight                   | 0,5 kg         |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,50 kg        |
| Gross weight approx.         | 0,50 kg        |
| Shipping weight              | 0,5 kg         |