

KERN 357-11

Individual weight, OIML Class M2, 1 kg, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-641
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	357

Measuring System

Test weight - nominal value	1 kg
Test weight - volume (cm ³)	119,0476
Test weight - volume uncertainty (cm ³)	2,4093
Permissible tolerance according to OIML (+/-)	160 mg
Test weight - Density [Min] (OIML)	2300
Weight OIML class	M2

Construction

Shape weight	Knob
Dimensions (Ø×H)	48×80 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	61×61×94 mm
Net weight	1 kg
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg