

KERN 357-12

Individual weight, OIML Class M2, 2 kg, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-642
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	357

Measuring System

Test weight - nominal value	2 kg
Test weight - volume (cm ³)	238,0952
Test weight - volume uncertainty (cm ³)	4,8186
Permissible tolerance according to OIML (+/-)	300 mg
Test weight - Density [Min] (OIML)	2300
Weight OIML class	M2

Construction

Shape weight	Knob
Dimensions (Ø×H)	60×100 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	72×72×115 mm
Net weight	2 kg
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg