

# KERN 357-13

Individual weight, OIML Class M2, 5 kg, Knob, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration	962-643
--------------------------------------	---------

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	357

## Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm <sup>3</sup> )	595,2381
Test weight - volume uncertainty (cm <sup>3</sup> )	12,0465
Permissible tolerance according to OIML (+/-)	800 mg
Test weight - Density [Min] (OIML)	2300
Weight OIML class	M2

## Construction

Shape weight	Knob
Dimensions (Ø×H)	80×147 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	100×100×165 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg