

KERN 357-14

KERN

Individual weight, OIML Class M2, 10 kg, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-644

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	357

Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm ³)	1190,4762
Test weight - volume uncertainty (cm ³)	24,093
Permissible tolerance according to OIML (+/-)	1,6 g
Test weight - Density [Min] (OIML)	2300
Weight OIML class	M2

Construction

Shape weight	Knob
Dimensions (Ø×H)	100×185 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	115×115×198 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg