

KERN 366-78

KERN

Individual weight, OIML Class M3, 20 kg, Block Eco, cast iron lacquered



Services (optional)

Article number for conformity assessment	952-645
Article number for DAkkS calibration	962-645

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

Measuring System

Test weight - nominal value	20 kg
Test weight - volume (cm ³)	2816,9014
Test weight - volume uncertainty (cm ³)	238,048
Permissible tolerance according to OIML (+/-)	10 g
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

Construction

Shape weight	Block Eco
Dimensions weight (W×D×H)	230×115×139 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×115×139 mm
Net weight	20 kg
Shipping method	Parcel service
Net weight approx.	20 kg
Gross weight approx.	20 kg
Shipping weight	20 kg