

KERN 367-05

Individual weight, OIML Class M3, 20 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-635
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	367

Measuring System

Test weight - nominal value	20 g
Test weight - volume (cm ³)	2,381
Test weight - volume uncertainty (cm ³)	0,0482
Permissible tolerance according to OIML (+/-)	25 mg
Weight OIML class	M3

Construction

Shape weight	Knob
Dimensions (Ø×H)	13×21,5 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	13×13×21 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg