

KERN 367-06

Individual weight, OIML Class M3, 50 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-636
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	367

Measuring System

Test weight - nominal value	50 g
Test weight - volume (cm ³)	5,9524
Test weight - volume uncertainty (cm ³)	0,1205
Permissible tolerance according to OIML (+/-)	30 mg
Weight OIML class	M3

Construction

Shape weight	Knob
Dimensions (Ø×H)	18×28 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	18×18×28 mm
Net weight	0,05 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg