

KERN 367-08

Individual weight, OIML Class M3, 200 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-638
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	367

Measuring System

Test weight - nominal value	200 g
Test weight - volume (cm ³)	23,8095
Test weight - volume uncertainty (cm ³)	0,4819
Permissible tolerance according to OIML (+/-)	100 mg
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

Construction

Shape weight	Knob
Dimensions (Ø×H)	28×47 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	28×28×47 mm
Net weight	0,2 kg
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg