

KERN 367-09

KERN

Individual weight, OIML Class M3, 500 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-639

Pictograms

STANDARD



1 DAY

Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	367

Measuring System

Test weight - nominal value	500 g
Test weight - volume (cm ³)	59,5238
Test weight - volume uncertainty (cm ³)	1,2046
Permissible tolerance according to OIML (+/-)	250 mg
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

Construction

Shape weight	Knob
Dimensions (Ø×H)	38×62 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	38×38×62 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg