

KERN ABP-A01

Ioniser to neutralise electrostatic charge

KERN



Category

| | |
|------------------|---------|
| Brand | KERN |
| Product category | Balance |
| Product group | Ionizer |

Construction

| | |
|-------------------|--------------|
| Dimension (W×D×H) | 124×87×49 mm |
| Material | plastic |

Environmental conditions

| | |
|---------------------------|--------|
| Storage temperature [Min] | -10 °C |
| Storage temperature [Max] | 40 °C |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 325×263×158 mm |
| Net weight | 0,85 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,85 kg |
| Gross weight approx. | 1,00 kg |
| Shipping weight | 2,7 kg |

Pictograms

STANDARD

