

# KERN MWA-A04

**KERN**

Handrail kit with folding patient seat for KERN MWA wheelchair platform scale 300K-1M



## Category

Brand	KERN
Product category	Balance
Product group	Tripod

## Construction

Dimension (W×D×H)	470×310×460 mm
Dimensions completely mounted (W×D×H)	1015×970×750 mm
Material	steel, lacquered
Stand height	1125 mm
Stand rotatable	✓

## Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	1000×1000×1816 mm
Net weight	8 kg
Shipping method	Freight
Net weight approx.	8 kg
Gross weight approx.	9 kg
Shipping weight	363,2 kg

## Pictograms

STANDARD

